

---

Subject: Re: Multiple DICOM images

Posted by [David Fanning](#) on Sat, 02 Feb 2002 00:14:37 GMT

[View Forum Message](#) <> [Reply to Message](#)

---

Akhila (idlfreak@yahoo.com) writes:

> Can anybody tell me how to read multiple dicom images - these are  
> continuous slices and are segmented. These images do not have any  
> patient details. Its just the image details and its an array[512,  
> 512].

Set the IMAGE\_INDEX keyword on READ\_DICOM to the index number of the image you are interested in reading. For example, if you want the fourth image in the stack:

```
image4 = Read_Dicom('mydicom.dcm', Image_Index=3)
```

Cheers,

David

--

David W. Fanning, Ph.D.

Fanning Software Consulting

Phone: 970-221-0438, E-mail: [david@dfanning.com](mailto:david@dfanning.com)

Coyote's Guide to IDL Programming: <http://www.dfanning.com/>

Toll-Free IDL Book Orders: 1-888-461-0155

---