
Subject: Re: Siemens Axiom Artis DICOM images
Posted by [Rick Towler](#) on Mon, 05 Jan 2004 19:50:30 GMT
[View Forum Message](#) <> [Reply to Message](#)

"Ed" wrote in message...

- > Dear group
- > Has anyone had any experience reading DICOM images taken from a
- > Siemens Axiom Artis Cardiac Cath Lab system?
- >
- > I am unable to get IDL to open any of the files we have written to CD
- > from this unit.
- >
- > READ_DICOM reports that:
- > The file sample.dcm is not in a supported DICOM format
- >
- > And when I attempt to access the image data using
- > array=obj->GetValue('7fe0'x, '0010'x, /NO_COPY)
- > array becomes -1
- >
- > I can access the data from a GE CT scanner, and from a GE Full Field
- > Mammo unit.
- >
- > I can also view all the images including the siemens ones using
- > Osiris.
- >
- > Any help would be much appreciated.

I suggest reading the IDL help on "DICOM, conformance summary".
Specifically the "Encapsulated Transfer Syntaxes Not Supported" section.

You will find that the Siemens files are in a format that IDL doesn't
support. I would suggest trying to export the images in a supported
transfer syntax (I don't know if Osiris can do this).

-Rick
