

---

Subject: Re: CT image without scaling

Posted by [David Fanning](#) on Thu, 24 Sep 2015 15:41:51 GMT

[View Forum Message](#) <> [Reply to Message](#)

---

Dete van Eeden writes:

> How can I display a CT image on idl without it being scaled to 0,255?

I'd try:

TV, image

Or, if the data is not in the range 0 to 255 already, scale it yourself with `BytScl` to whatever range you want within the values 0 to 255 required to put an image on the display.

Cheers,

David

--

David Fanning, Ph.D.

Fanning Software Consulting, Inc.

Coyote's Guide to IDL Programming: <http://www.idlcoyote.com/>

Sepore ma de ni thue. ("Perhaps thou speakest truth.")

---