
Subject: IDLffDICOMEx

Posted by [hyperjiver](#) on Wed, 31 May 2006 21:27:12 GMT

[View Forum Message](#) <> [Reply to Message](#)

Hi..

I tried using IDLffDICOMEx command to clone a file.

The syntax I used is given here:

```
testobj = obj_new('IDLffDICOMEx', fname, CLONE = tempdcmfile[0])
```

I get this error message....

```
% OBJ_NEW: Error: Failed to reset filename when cloning existing file
, (null), Null parameter
% Execution halted at: FMRIIMAGE_DICOM_CONVERT_GRAY 240
D:\RSI\CCHIPS.2005\devlib\display_utils\fmriimage\fmriimage_dicom_convert_gray.pro
%          FMRIIMAGE_EVENT 331
D:\RSI\CCHIPS.2005\devlib\display_utils\fmriimage\fmriimage_perfusion.pro
%          WIDGET_PROCESS_EVENTS
%          $MAIN$
```

can anybody help??

Thanks,
HJ
